


Cal Poly Pomona, L.S. STARKEN NOTEBOOK SAMPLE

Title: Synthesis of Cyclohexene 1/16/17

Purpose: The dehydration of cyclohexanol w/ H₂SO₄ will produce cyclohexene.



Reagent*	MW	Dens.	Am't. (2.1 mL)	Mol	Equiv.*	Remarks*
cyclohexanol	100.16	0.96	2.0g	0.020	1.0	bp 161°C
H ₂ SO ₄ (conc.)	98.07	1.84M	0.4 mL	0.007	0.37	caution! Corrosive
water	18.0	1.0	0.3 mL	-	-	-

cyclohexene: 82.14 0.81 1.64g* 0.020 1.0 bp 83°C, flammable!

Workup materials: toluene, ag. NaCl soln, CaCl₂

Calculations*

$2.0g \text{ cyclohexanol} \times \frac{\text{mol}}{100.16g} = 0.020 \text{ mol}$

theor. yd: $0.020 \text{ mol cyclohexene} \times \frac{82.14g}{1 \text{ mol}} = 1.64g$

Product is final entry in table

Sketch of Apparatus*




FIG 1 Fractional Distillation Apparatus

Workup Outline

Step	Purpose
A. Wash with ag. NaCl	remove H ₂ O (dry)
B. Add CaCl ₂	remove trace H ₂ O (dry)
F. Rotovap	remove ether

page 1...

Procedure*

- 1) To a 5 mL R.B flask, add the H₂O, then the H₂SO₄ and then the cyclohexanol + boiling chip
- 2) Cover (?) + shake to mix layers. (heat?)
- 3) Assemble apparatus in Fig 1.
- etc

Workup:

- 1) Transfer distillate to a test tube
- 2) Wash w/ equal vol. of satd NaCl soln.
- 3) Remove aq. layer (bottom)
- 4) Dry org. layer (top) w/ CaCl₂ (product where!)

Waste Disposal

Aqueous solns flush down drain...

Observations
 → alcohol is clear, colorless liq.
 → Turned bright pink.
 → flask got hot.
 → distilled @ 80°C
 ...
 etc.
 → spilled some here
 → used Na₂SO₄ as the drying agent
 → obtained 0.02g product as a clear, colorless liquid
 *to be written in ink during lab!

*Notes:

Prelab is worth 50% of total lab grade!
 If it is not done before lab, 25 pts are deducted.

Prelab Assignment

Reagents list all species used in rxn (include solvents) but not workup materials.
 Equiv. is molar ratio or # of equivalents used.
 Remarks should include mp (solids), bp (liquids), solubility, + any relevant hazards/precautions
 1.64g This is the theoretical yield of product, as calculated.
 Calculations should include proper units
 Sketch new apparatus only if used before, reference that sketch's page # here. (ie, see Sep. Funnel sketch on p15)
 Procedure Brief, step-by-step outline. Indicate where the reaction workup (product isolation) begins. Finish with notes on waste disposal. Modify as needed while conducting experiment!

page 2...